

Overview

- SPM Components
- Spatial Transformations
- Model specification
- Comments

1

SPM: Overview

- Library of Matlab and C functions
- “Point-&-click” graphical user interface
- Consists of 4 general components
 - Preprocessing
 - Model Specification & Fitting
 - Inference & Results Interrogation
 - Supplemental Tools

2

SPM: Preprocessing

- Elimination of systematic variation before statistical modeling
- “Realign”ment
 - Intrasubject registration
 - Motion correction
- “Coregister”ation
 - Intrasubject, *inter*modality registration
 - Registration of MR images with different TR/TE
 - PET-MR registration
- Spatial “Normalize”ation
 - *Intersubject* registration
 - Register subject anatomy to atlas space

3

SPM: Preprocessing

- Spatial “Smooth”ing
 - Blur data with Gaussian filter
 - Smoothing data into submission...
 - To satisfy random field theory assumptions
 - For intersubject analyses
- “Segment”ation into GM/WM/CSF
 - Usually not directly used
 - Useful for structural studies

4

SPM: Model Specification & Fitting

- “fMRI models”
 - Specify a fMRI design
 - Review a fMRI design
 - Estimate a fMRI design
- “Explore design”
 - Examine correlation of predictors
 - Power spectrum of experimental effects
- “Basic models”
 - t-tests (Two sample, One sample, Paired)
 - Regression

5

SPM: Inference & Results Interrogation

- “Results” button
- First brings up “Contrast Manager”
 - Can define single (t) or sets (F) of contrasts
- Then shows MIP in Graphics window
 - MIP = Maximum Intensity Projection
 - Glass Brain
 - Can “surf” by dragging cursor

6

SPM: Inference & Results Interrogation

- Interactive window
 - p-values
 - Corrected for whole brain or subregion
 - Plotting of time courses
 - “Overlays”
 - Superimpose results on other images
 - Current location and value

7

SPM: Miscellaneous Tools

- “Display”
 - Display an image with orthogonal sections.
 - Check intensity values
 - Change origin
 - Change world space
 - I.e. apply rotations/translations to mat file
- “Check Reg”
 - Display multiple images
 - *Essential* tool for assessing estimated or assumed alignment of images.
 - All images are displayed in the space of the first image.

8

SPM: Miscellaneous Tools

- “ImCalc”
 - Image calculator
 - Give one or more images, perform Matlab arithmetic and write out result
- “HDR edit”
 - Edit Analyze headers
 - Can change just a few fields, leaving others untouched.
 - Can apply changes to multiple images.
- “Utils”
 - Change directory
 - Results are written to current directory!
 - delete files, etc.

9

SPM: Miscellaneous Tools

- “Rendering”
 - Sexy brains
- “Defaults”
 - Change defaults for one session.
 - To make permanent changes, edit `spm_defaults.m`
- “Help”
 - Some help
 - Also not “Tool tip” help balloons

10

SPM: Pipeline

-

11

Spatial Transformations: Image File

- Data must be in analyze format
- `.img` Raw, binary data; single 3D volume
- `.hdr` Short binary header
 - Image dimension
 - Voxel size
 - Origin, in voxels
 - First element 1, not 0
- `.mat` Optional, SPM extension
 - Defines transformation from voxel to world space
 - If exists, `.hdr` voxel size & origin are ignored
 - Origin can be represented as mm location
 - e.g. between voxels

12

Spatial Transformations: Voxel vs. World Space

- Voxel Space
 - Just the original image
 - No reorientations or flips
- World Space
 - Space defined by `.mat` file
 - `.mat` - 4×4 homogeneous transformation matrix M
 - Let v be a voxel location indexed from (1,1,1)
 - Then $w=M*[v;1]$ is that location in world space
 - Has units of *mm*
 - Can represent rotations, translations and flips

13

Spatial Transformations: Coregistration & Realignment

- Coregistration & Realignment are rigid body transformations
 - Subject's doesn't change size or warp between scans
 - Well, actually...
- Each requires a "Target" and an "Object"
 - Target — Fixed image
 - Object — Image that is transformed
- SPM modifies the `.mat` file of the *object* image
 - Unless you explicitly ask it to, it doesn't write out an image

14

Spatial Transformations: Coregistration & Realignment

- This is very cool!
 - You have 1000 images
 - You realign 2nd through the 1000th to the 1st
 - Target=1st, Object=2nd-1000th
 - Instead of doubling the disk space required by writing out 1000 images, you only need to modify the `.mat`'s!

15

Spatial Transformations: Coregistration

- `.mat`'s handy for coregistration
- e.g. Coregister
 - Target: Functional space, 7 slices
 - e.g. 7 slice, low resolution anatomical in space of func's, "t1_gre"
 - Object: Structural space
 - e.g. high resolution anatomical, 120 slices difference space, "t1_spgr"
- Reslicing an image, creating "t1_spgr-gre"
 - You get a 7 slice, low resolution SPGR

16

Spatial Transformations: Coregistration

- But don't have to reslice
 - `t1_spgr`'s world space now corresponds to `t1_gre`
 - The `.mat` file has captured the transformation *without* throwing away any information

17

Spatial Transformations: Easy to forget about `.mat`'s

- Easy to forget about `.mat`'s
 - If you want to start over and re-do registration *must delete the .mat's!*

18

Spatial Transformations: Left is Right is Left is Right is ...

- Two conventions for images
- Neurological
 - On the screen, Left is Left side of subject
 - As if standing behind the patient
- Radiological
 - On the screen, Left is Right side of subject
 - As if standing at the foot of the patient
- Standard in clinical radiology is, um, radiological
- SPM uses Neurological
 - Be very careful!

19

Spatial Transformations: Left is Right is Left is Right is ...

- Best to make a lab convention
 - fMRI Lab will (*will*) produce everything in neurological
- If paranoid
 - Tape a vitamin E capsule to head and call me in the morning

20

SPM: Specifying fMRI Models

- fMRI design specification
 - Most nonintuitive thing about SPM (IMHO)
- SOA
 - Stimulus Onset Asynchrony
 - Spacing between *beginning* of events
 - A “Fixed” or regular spacing means only need to specify the onset time of first event/block
 - A “Variable” SOA means you’ll need to specify a vector of onset times

21

SPM: Specifying fMRI Models

- Onsets
 - Offset time of each block/event start
 - Units of *scans*
 - Offset of 0 = first scan
- Events vs Blocks
 - Just a difference in stimulus train
 - Events - Stick functions
 - Blocks - Box Cars
 - Or can specify a basis of stimulus trains
 - e.g. Half sinewaves, or Sine/Cosine basis

22

SPM: Specifying fMRI Models

- “Volterra” Interaction among trials
 - An attempt to model nonlinear interactions between trials with a linear basis
 - Not for the weak of heart

23

SPM: SPMese Decoder Ring

- “Spatial Normalization”
 - Intersubject registration
- “Coregistration”
 - Intermodality registration (intrasubject)
- “Registration”
 - Intrasubject registration (movement corr)
- “Results”
 - Inference
- “Interscan Interval”
 - TR
- “SOA”
 - Stimulus Onset Asynchrony

24

SPM: Perspective

- SPM tries to be a “sole source”, turn key solution for fMRI
- At this point, *there can be no such thing!*
 - This is impossible in a cutting edge, rapidly evolving field where each dataset has 10's of millions of observations!!!
- Don't let SPM be a black box!
- Understand what each component does
- Understand how to get at the data
 - e.g. using 'Display' and 'Check Reg' alot.